KIRWAN SURGICAL PRODUCTS LLC
INSTRUCTIONS FOR USE –
DISPOSABLE BIPOLAR CORDS

This product is disposable, supplied sterile and intended for single use only. Prior to the expiration date, contents are sterile if the package has not been opened or damaged. Do not resterilize.

INTENDED USE: Disposable cords are designed to connect an electrosurgical device to an electrosurgical generator. Please refer to the labeling to determine the particular connectors and specific use of the cord.

CONTRAINDICATIONS: Any use of this instrument for tasks other than for which it is indicated will usually result in a damaged or broken instrument.

SETUP AND USE

WARNING: CONNECT BIPOLAR ACCESSORIES TO THE BIPOLAR RECEPTACLE AND MONOPOLAR ACCESSORIES TO THE MONOPOLAR RECEPTACLE ONLY. IMPROPER CONNECTION OF ACCESSORIES MAY RESULT IN INADVERTENT ACCESSORY ACTIVATION OR OTHER POTENTIALLY HAZARDOUS CONDITIONS.

WARNING: DO NOT REUSE OR RE-STERILIZE. PATIENT CROSS CONTAMINATION AND/OR PRODUCT MALFUNCTION MAY OCCUR

Attach the sterile cord to the sterile instrument, ensuring that the contact pins are fully seated in the cord receptacles. Connect the cord to the generator only while it is in the "OFF" or "STANDBY" mode. Failure to do so may result in injury or electrical shock to the patient or operating room personnel.

At the lowest power setting, test the instrument connections by pressing the generator's activating switch. If the generator fails to activate, check the connection with the cord.


Power setting guidelines may vary due to differences in surgical techniques, patients, electrodes and surgical set-up. Start at the lowest power setting and increase as necessary to achieve the desired clinical effect.

CAUTION: BECAUSE OF THE VARIABILITY OF OUTPUT VOLTAGES AND MODES FROM GENERATOR TO GENERATOR, DO NOT USE THE BIPOLAR CORD WITH GENERATOR SETTINGS HAVING A BIPOLAR OUTPUT VOLTAGE EXCEEDING 1000 Vp-p. REFER TO THE APPROPRIATE ELECTROSURGICAL GENERATOR MANUAL FOR INDICATIONS AND INSTRUCTIONS ON VOLTAGE OUTPUT CHARACTERISTICS TO ENSURE THAT ALL SAFETY PRECAUTIONS ARE FOLLOWED.

	SEE INSTRUCTIONS FOR USE
	USE BY: YEAR-MONTH
	SINGLE USE – DO NOT REUSE
	Do not resterilize
	BATCH CODE
	CATALOG NUMBER
	STERILE – STERILIZED BY IRRADIATION

Sterility guaranteed unless package opened or damaged. Do not resterilize.

NO LATEX Not made with natural rubber latex.

 CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

CE 1639

Kirwan Surgical Products LLC
 180 Enterprise Drive
 Marshfield, MA 02050, USA
 Tel: 888-547-9267
 Product of USA

Medical Product Service
 Borggasse 20
 35619 Braunfels, GERMANY

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