


Disposable Bipolar Suction Coagulators

| ITEM # | TOTAL LENGTH | DESCRIPTION | FEATURES |
|---------|-----------------|------------------------------------|----------|
| 14-9000 | 7 1/2" (19.1cm) | 8 French | IN |
| 14-9001 | 7 1/2" (19.1cm) | 11 French | IN |
| 14-9002 | 5 1/2" (14cm) | 13 French | IN |
| 14-9004 | 5 1/2" (14cm) | 15 French | IN |
| 14-9006 | 7 1/2" (19.1cm) | 13 French | IN |
| 14-9008 | 7 1/2" (19.1cm) | 15 French | IN |
| 14-9020 | 5 1/2" (14cm) | 15 1/2 French, Angled/Chiseled Tip | IN |
| 14-9030 | 11" (27.9cm) | 8 French | IN |
| 14-9031 | 11" (27.9cm) | 11 French | IN |
| 14-9032 | 11" (27.9cm) | 13 French | IN |
| 14-9034 | 11" (27.9cm) | 15 French | IN |
| 14-9041 | 7 1/2" (19.1cm) | 8 French, Angled Up | IN |
| 14-9042 | 7 1/2" (19.1cm) | 8 French, Angled Down | IN |
| 14-9043 | 7 1/2" (19.1cm) | 11 French, Angled Up | IN |
| 14-9044 | 7 1/2" (19.1cm) | 11 French, Angled Down | IN |



Disposable Monopolar Suction Coagulators with 10' Resectoscope Pin Cord

| ITEM # | TOTAL LENGTH | DESCRIPTION | FEATURES |
|---------|---------------|-------------|----------|
| 30-5200 | 5 1/2" (14cm) | 10 French | IN |
| 30-5201 | 8" (20.3cm) | 10 French | IN |